# 《Introduction to tissue preparation instruments》

## Introduction

Advantages of each frozen sections and paraffin sections.

Morphology → Paraffin section  
Proteins, Lipids, Preservation of antigenicity → Frozen section

## Location

Cryostat :

CM1900	room206
CM1950	room511

Equipments for paraffin specimens :

HISTOS5, EG1160, RM2125RT	room511
RM2225	room206
AS-200H	Liaison Laboratory

## 【Tips for specimen preparation】

### <Frozen specimens>

Fixed or unfixed tissue → Frozen  
→ Sectioning → Staining → Microscope

CM1950


CM1900, CM1950 :  
Cryostat.

When PBS and fixative solution remains in a block, slicing is difficult. Please rinse with compound once and transfer samples to a new compound for making blocks.

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## <Paraffin specimens >

Fixed tissue → Paraffin Processing  
→ Sectioning → Staining → Microscope

1.HISTOS5

**HISTOS5** : Microwave Histoprocessor.  
(dehydration, delipidation ~ replacement with paraffin).

It is necessary to fix strongly for a microwave processing.

Fixation: 10% formalin, room temperature, 3 days

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Please confirm clamping levers in three places are tightened firmly, or slices could be alternate or uneven thickness.

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AS-200H

### AS-200H : Auto Slide Preparation System

This is useful when preparing a lot of specimens from a relatively large sample. (example:slicing whole mouse embryos and adult organs)

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